Shown to heal 2-3x faster than Standard of Care¹

NEXT SCIENCE°

30 mL 1 oz.

NEXT SCIENCE

BLASTX®

Antimicrobial Wound Gel

BLASTX Antimicrobial Wound Gel aids in wound management by maintaining a moist wound environment. While in place, the antimicrobial properties of the gel inhibit the growth of microorganisms in the product.

LEARN MORE AT NEXTSCIENCE.COM/BLASTX





linical Benefits

- Provides a moist wound environment conducive to healing²
- Biocompatible³
- Prevents bacterial growth within the gel4
- Broad spectrum efficacy⁵



Indication for BLASTX

BLASTX is indicated for the management of wounds such as stage I-IV pressure ulcers, partial- and full-thickness wounds, diabetic foot and leg ulcers, post-surgical wounds, first and second degree burns, and grafted and donor sites.




Instructions for use

- 1. Apply BLASTX directly to the wound at a recommended thickness of 3mm.
- 2. Cover with an appropriate dressing.
- 3. After opening, product may be used until labeled expiration date.



Successful Wound Management

BLASTX Antimicrobial Wound Gel Demonstrates Successful Management of the Wound²



A successfully treated wound is defined as >50% reduction of wound volume in 4 weeks. Standard of Care (SOC) is custom antimicrobial cream.

*p-value = 0.05




🔙 Case Study

Patient had a history of diabetes mellitus and peripheral artery disease. Before using BLASTX, patient received negative pressure treatment for two months and healing stalled. BLASTX started on April 23rd and healed in 16 days.

- 63-year-old male
- Necrotizing cellulitis and deep abscess





Wound healed after 16 days

- 1. Kim D, Namen II W, Moore J, Buchanan M, Hayes V, Myntti MF, Hakaim A. Clinical Assessment of a Biofilm-disrupting Agent for the Management of Chronic Wounds Compared With Standard of Care: A Therapeutic Approach. Wounds 2018 May;30(5):120-130. Epub 2018 Jan 26. PMID: 29521641.
- 2. Wolcott, R. (2015). Disrupting the biofilm matrix improves wound healing outcomes. Journal of Wound Care. 24(8),366-71. doi: 10.12968/jowc.2015.24.8.366
- Satisfies ISO 10993-1
- 4. As ner USP 51
- 5. Data on file: TR-03-16-006

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CONTAINS Benzalkonium Chloride 0.13%, Polyethylene Glycol 400, Polyethylene Glycol 3350, Sodium Citrate, Citric Acid, and Water.

CONTRAINDICATIONS BLASTX Antimicrobial Wound Gel should not be used if there is a history of allergy to any of the ingredients.

PRECAUTIONS For external topical use only. Do not ingest. Must not be used in or near the eyes. Do not cover with alginate dressings



